# PREP KIT

INSTRUCTION FOR USE KULLANMA KILAVUZU

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- 1 Unit of Needle (30G – 13 mm)
- 1 Unit of Needle (21G – 40 mm)
- 1 unit of Blood Collection Needle

#### Intended Use of the Product

The device is designed to be used for the safe and rapid preparation of autologous platelet-rich-plasma (PRP) from a small sample of blood at the patient's point of care. The PRP is prepared prior to application to a local area of soft tissue or hard tissue (including bone) site as deemed necessary by the clinical use requirements. And its intended use is for the application of PRP in order to achieve the stimulation and acceleration of soft tissue and hard tissue (including bone).

#### Intended Patient Population:

is an autologous application technique therefore patient applications are inherently safe from infectious diseases such as HIV, Hepatitis, Creutzfeldt-Jakob disease.

This application can be used safely and effectively on all patients. The safety and effectiveness of PRP treatment have not been established in children and pregnant women.

The device is intended to be used for the treatment of the following conditions:

- alopecia areata,
- joint pain (particularly in patients with osteoarthritis of the knee),
- rotator cuff injuries,
- and lateral epicondylitis (tennis elbow).

#### The Performance of the Device Intended by the Manufacturer

PRP (Platelet Rich Plasma) is used for many years in the frame of regenerative medicine, which allows the cellular healing on bone and soft tissue recovery. Bone and soft tissue healing is occurred by the platelet concentration and re-injection locally to the related part of the body.

Currently accepted technique for avoiding coagulation of the whole blood, referred to Sodium Citrate. The procedure uses anticoagulated blood, collected by routine phlebotomy. Buffered Citrate Solution blocks the coagulation by mixing gently the whole blood together with 3,20% concentrated 0,1M of Sodium Citrate. The prepared material is then centrifuged to obtain PRP. The recommended G Force for concentration is given in the "Procedure" section according to the centrifuge rotor. The blood separation system takes the advantage of centrifugation force effecting on the weight of the cells. The separation occurs during centrifugation when the denser blood components move due to a gravity force. The PRP Tube shall not be mixed after the centrifugation. According to the centrifuge rotor used and the desired protocol, PRP is harvested carefully above the erythrocytes (Red Blood Cell layer). The collected PRP is moved to Re-suspension tube for re-suspension process.

The result is a convenient tube system, collection of the whole blood and

concentration of Platelet Rich Plasma. PRP Tubes included in PRP KIT are interior sterilized using irradiation. PRP KIT reduces the risk of air contamination, eliminates the need of any other tubes or pipettes.

#### Clinical Safety Information

The safety and efficacy of this device has been proven by Risk Analysis Report and Clinical Evaluation Report, attached to the Technical File of PRP Kit.

#### Residual Risk and Side Effects Information

According to the risk assessment, risk management, risk report procedures of this device, it shows no residual risks. With respect to the GSPR 3, the risk management report showed that there were no serious side effects and the performance of the product was verified by the clinical evaluation report.

#### Intended User

This product should be used by physicians or physician directed allied health care professionals with adequate training in the use of the device.

#### Procedure

Apply the following instructions to obtain Platelet Rich Plasma by using PRP KIT:

- Open the box.
- Open the sterile blister package.
- Attach Blood Collection Needle to Holder.
- Decontaminate the skin and then, using a needle to draw blood, insert into the appropriate vein.
- Attach PRP Tube inside the holder to Blood Collection Needle. (Routine Phlebotomy)
- PRP tube is a vacuumed tube. The tube collects enough amount of blood and stops itself. Use both tubes without removing the Blood Collection Needle from the vessel, in order to obtain the optimal amount of autologous material for the procedure.
- After blood collection, place both tubes in a centrifuge and centrifuge the material as recommended in the table below. The choice of procedure depends on the individual indications of the patient and the recommendations of the professional.

	Fixed rotor centrifuge	Swing rotor centrifuge
Pure PRP Protocol	500 RCF & 3 Minutes	1150 RCF & 4 Minutes
Amber PRP Protocol	500 RCF & 6 Minutes	1150 RCF & 8 Minutes
	500 RCF=2000 RPM, 1150 RCF=2900 RPM	
	AMBER PRP: LR-PRP (Leukocyte Rich PRP)	
	PURE PRP: LP-PRP (Leukocyte Poor PRP)	

- Buckets shall be compatible for PRP Tubes (16 X 100 mm tube size).

## Instruction For Use

EN

This Instructions for Use document covers all brands of the product with reference 1000.

#### The Name or Trade Name of the Medical Device

PRP KIT

The product is a PRP KIT, sterile, single use medical device.

#### The Name and the Address of the Manufacturer

T-Biyoteknoloji Laboratuvar Estetik Medikal Kozmetik San. ve Tic. Ltd. Şti. Dumlupınar Mah., Kuleli Cad., No:3 Nilüfer – Bursa / TÜRKİYE  
Manufacturer ID No (SRN): TR-MF-000016403  
Phone: +90 224 246 85 22  
E-mail: info@tlab.com.tr

#### Product Content List

- 2 units of PRP Tube(s)
- 1 unit of Re-suspension Tube
- 1 unit of Long Needle (18G– 88mm)
- 1 unit of Holder
- 1 unit of 1ml Injector
- 1 unit of 5ml Injector

- Attach Long Needle to 5 mL Luer-Lok Syringe.
- Collect the PRP by sticking the Long needle into the PRP Tube above the RBC Layer (Do not open the Lid of the tube unless it is performed under Laminar Flow Conditions)
- After the Pure PRP Protocol, collect all plasma above the erythrocyte layer.
- After the Amber PRP Protocol, collect the buffy coat along with 2ml plasma above the erythrocyte layer.
- Collect carefully PRP Layer from first PRP Tube.
- Carry PRP to the Re-suspension Tube.
- Apply again the same protocol for the other PRP Tube.
- Carry again PRP to the Re-suspension tube.
- Gently shake the Re-suspension tube for 30 seconds to homogenize the settlement of the possible clumped cells in the plasma milieu. (Re-suspension process)
- Attach the application needle (Green Needle - 21G) to 1ml Luer-Lok syringe.
- Collect the re-suspended PRP inside from the Re-Suspension Tube.
- Change application needle (green needle - 21G) to injection needle (yellow needle - 30G)
- Use an appropriate injection method. From the time of centrifugation, the plasma should be used as soon as possible.

#### Storage, Shelf Life and Handling Information

- Do not expose to direct sunlight.
- Do not use if the inner sterile blister package is damaged or opened.
- Expiration date is available on the label.
- Store at temperatures between 5°C and 30°C.

#### Precautions, Warnings, Measures to Be Taken and Limitations

- Use proper safety precautions to guard against needle sticks.
- Follow manufacturer instructions when using centrifuge.
- Do not use sterile component of this kit if package is opened or damaged.
- Do not open the lids of the tubes.
- Single use device, do not reuse.
- The user should be thoroughly familiar with the equipment and the surgical procedure prior to using this device.
- The patient should be made aware of general risks associated with treatment and possible adverse effects.
- Do not dispose of this product in the unsorted municipal waste stream.
- Long Needle is not intended to be used on human. It is just intended to be used to collect PRP Layer inside from the PRP Tubes).

#### Possible Adverse Events

- Damage to blood vessels, hematoma / delayed wound healing and/or infection associated with blood draw.
- Temporary or permanent nerve damage that may result in pain or numbness.
- Early or late postoperative infection.

#### Possible Contraindicative Situations

Platelet Dysfunction Syndrome, Critical

Thrombocytopenia, Hypofibrinogenemia, Hemodynamic instability, Auto-immune diseases, Malignancy, Sepsis, Acute Chronic infections, Chronic Liver Pathology, Anti-coagulation Therapy (ASA therapy etc.), Pregnancy, Newborns, Not recommended for use under 18 years of age.

#### Symbols and Explanations

	CE Number Notified Body: Suzutec (2195)		Do Not Re-use
	Manufacturer		Do Not Re-sterilize
	Country of Manufacture & Date of Manufacture		Do Not Use if Package is Damaged
	Use by date		Keep Away from Sunlight
	Catalogue Number		Keep Dry
	Batch Code		Temperature Limit
	Sterilized Using Irradiation		Consult Instructions for Use or Consult Electronic Instructions for Use
	Sterilized Using Ethylene Oxide		Caution
	Single Sterile Barrier System		Unique Device Identifier
	Medical Device		Non-Pyrogenic
	Importer		Distributor

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