TR
EN
EL
DE
DK
FR
ITA
SK
ESP

PRP

SYRINGE

INSTRUCTION FOR USE

TR	2
EN	4
EL	8
DE	11
DK	15
FR	18
ITA	21
SK	24
ESP	28
SYMBOLS	32

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Instruction for Use EN

This Instructions for Use document covers all brands of the product with reference 500.

THE NAME OR TRADE NAME OF THE MEDICAL DEVICE

PRP SYRINGE

The product is a sterile, single use medical device.

THE NAME AND THE ADDRESS OF THE MANUFACTURER

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E-mail: info@tlab.com.tr

PRODUCT CONTENT LIST

- 1 unit of 10ml PRP SYRINGE
- 1 unit of Luer Female Cap
- 1 unit of Luer-Luer Female Adapter

INTENDED USE OF THE PRODUCT

The device is designed to be used for the safe and rapid preparation of autologous platelet-rich-plasma (PRP) from a small sample of blood at the patient's point of care. The PRP is prepared prior to application to a local area of soft tissue or hard tissue (including bone) site as deemed necessary by the clinical use requirements. And its intended use is for the application of PRP in order to achieve the stimulation of soft tissue and hard tissue (including bone) and acceleration of soft tissue and hard tissue (including bone) healing.

INTENDED PATIENT POPULATION

Platelet rich plasma is an autologous application technique therefore patient applications are inherently safe from infectious diseases such as HIV, Hepatitis, Creutzfeldt-Jakob disease. This application can be used safely and effectively on all patients. The safety and effectiveness of PRP treatment have not been established in children and pregnant women.

The device is intended to be used for the treatment of the following conditions:

- alopecia areata,
- joint pain (particularly in patients with osteoarthritis of the knee),
- rotator cuff injuries,
- lateral epicondylitis (tennis elbow)
- and Chronic Plantar Fasciitis

THE PERFORMANCE OF THE DEVICE INTENDED BY THE MANUFACTURER

Platelet cells have been used for many years in the field of regenerative medicine for cellular healing of bone or soft tissue. Bone and soft tissue healing is achieved by re-injecting the obtained platelet concentration locally into the relevant area of the body.

Platelet concentration is very important for PRP/PRF

treatments. Blood drawn from the patient is centrifuged to concentrate the platelet layer. The recommended G-force for concentration is 1500 G for 4 minutes. The blood separation system makes use of the centrifugal force that affects the weight of the cells. Separation occurs during centrifugation when dense blood components move due to gravitational force.

PRP Syringe should not be mixed after centrifugation. A 3–4 ml layer is carefully transferred to another syringe down to the erythrocyte layer (Red Blood Cells).

As a result, the device is a suitable PRP/PRF syringe system for blood collection and Platelet-Rich Plasma concentration. PRP Syringe reduces the risks of ambient air contamination and eliminates the need for other accessories such as tubes or pipettes.

CLINICAL SAFETY INFORMATION

The safety and efficacy of this device has been proven by Risk Analysis Report and Clinical Evaluation Report, attached to the Technical Construction File of PRP SYRINGE.

RESIDUAL RISK AND SIDE EFFECTS INFORMATION

According to the risk assessment, risk management, risk report procedures of this device, it shows no residual risks. With respect to the requirements of Essential Requirements 6, the risk management report showed that there were no serious side effects and the performance of the product was verified by the clinical evaluation report.

INTENDED USER

This product should be used by physicians or physician-directed allied health care professionals with adequate training in the use of the device.

PROCEDURE

To obtain PRP/PRF using a syringe, follow these steps in order.

- Open the box.
- Open the sterile blister pack.
- Insert a hypodermic or butterfly blood collection needle into the syringe.
- Prepare the patient's arm in order to draw blood (sanitize the skin).
- Insert the needle into the cleaned area of the patient's arm.
- Prepare the balance weight to run the centrifuge smoothly.
- Collect blood up to the 10ml mark on the syringe. (Phlebotomy).
- Remove the blood collection needle from the syringe and close the syringe with the locked cap inside the blister pack.
- Remove the syringe piston rod by turning it and place it in the centrifuge without delay.

- Place the syringe in the centrifuge with the cap facing upward.
- Set the centrifuge to 1500G for 4 minutes and start the centrifuge.
- After centrifugation, remove the syringe from the centrifuge, do not shake or rotate it, and hold it upright with the cap on top.
- Reinsert the piston rod into the syringe.
- Open the syringe cap, and while holding the syringe upright, attach the locked connector to the syringe.
- Attach another injector to the locked connector.
- During the transfer of the prepared PRP/PRF layer from the syringe, press the piston until it reaches the blood cells.
- Do not collect red blood cells (RBCs).
- The final product is ready for injection.

STORAGE, SHELF LIFE AND HANDLING INFORMATION

- PRP SYRINGE has been sterilized with Ethylene Oxide.
- Do not expose to direct sunlight.
- Store in a cool, dry place.
- Do not use if the inner sterile blister package is damaged or opened.
- Product shelf life is 5 years. Use by Date is indicated on the product label.
- Store at temperatures between 5°C and 30°C.

PRECAUTIONS, WARNINGS, MEASURES TO BE TAKEN AND LIMITATIONS

- Use proper safety precautions to guard against needle sticks.
- Follow manufacturer instructions when using centrifuge.
- Do not use sterile component of this kit if package is opened or damaged.
- Single use device, do not reuse.
- The user should be thoroughly familiar with the equipment and the surgical procedure prior to using this device.
- The patient should be made aware of general risks associated with treatment and possible adverse effects.
- The safety and effectiveness of PRP treatment have not been established in children and pregnant women.
- Do not dispose of this product in the unsorted municipal waste stream.

POSSIBLE ADVERSE EVENTS

- Damage to blood vessels, hematoma / delayed wound healing and/or infection associated with blood draw.
- Temporary or permanent nerve damage that may result in pain or numbness.
- Early or late postoperative infection.
- Inflammatory reaction.

POSSIBLE CONTRAINDICATIVE SITUATIONS

Platelet Dysfunction Syndrome, Critical

Thrombocytopaenia, Hypofibrinogenaemia, Hemodynamic instability, Auto-immune diseases, Malignancy, Sepsis, Acute Chronic infections, Chronic Liver Pathology, Anti-coagulation Therapy (ASA therapy etc.), Pregnancy, Newborns

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